

EEG, MEG and neuromodulatory approaches to explore cognition: Current status and future directions.

Carolina Beppi, Inês Ribeiro Violante, Gregory Scott, Stefano Sandrone.

Presented By :
Deepa Tilwani
AI Institute
University of South Carolina
Nov, 22nd 2021

Link to paper : <https://www.sciencedirect.com/science/article/pii/S0278262620302803?via%3Dihub>


Paper Focus Points:

- How non-invasive brain stimulation (NIBS) methods might better apply to certain cognitive functions or understand their relative efficacy for improving cognitive function ?

Cognitive Function


Six principle domain of Cognitive Functions:

1. Complex attention
2. Executive Function
3. Learning and Memory
4. Language
5. Motor Function
6. Social Cognition



EEG : Introduction

- **Electroencephalogram (EEG) electrodes**
- Scalp recording of **electrical activity** of cortex => waveform signals
- **Microvolts (μV)** – small!
- Role of EEG in neuroimaging:
 - Identify **neural correlates**
 - **Diagnose** epilepsy, sleep disorders, anaesthesia, coma, brain death



MEG: Introduction



<http://www.admin.ox.ac.uk/estates/capitalprojects/previouscapitalprojects/megscanner/>

- **Magnetoencephalography**
- Direct external recordings of **magnetic fields** created by electrical currents in cortex
- Measured in **fT** to **pT**
- Role of MEG in neuroimaging:
 - **Neural correlates** of cognitive/perceptual processes
 - **Localise** affected regions before surgery determine regional and network functionality
- **Excellent spatial resolution**
good for functional mapping of specific cortex (M1, V1) during behavioural, cognitive, perceptive tasks
- **Surgical planning** in patients with brain tumours or intractable epilepsy
- **Research** into whole-brain network connectivity
Millisecond temporal resolution

EEG: Frequency Spectrum

- **5-50 μV** , mostly below 30 μV
- Sharp spike-waves, **light sleep** stages

Beta (β) 13-30 Hz

Frontally and parietally



- **5-120 μV** , mostly below 50 μV

- **Awake**, eyes closed, mental inactivity, physical relaxation

Alpha (α) 8-13 Hz

Occipitally



- **20-200 μV**

Theta (θ) 4-8 Hz

Children, sleeping adults



- **5-250 μV**

- Abnormality in waking adults, accompaniment of **deep sleep**

Delta (δ) 0.5-4 Hz

Infants, sleeping adults



+ **Gamma waves**



31-100 Hz, 10 μV

‘binding of consciousness’, unity of perception

Spikes 3 Hz

Epilepsy - petit mal


V [μV]
200
100
0



Major brain rhythms as classified by their frequency span.

Non-Invasive Brain Stimulation (NIBS) Techniques

- Different from the recording of brain (EEG,MEG,fMRI)
- NIBS are causations methods (which causes change in the brain) and then we measure changes in behaviour or cognition.
- Can help studying neuronal changes in ms.
- Able to detect how much effect on brain or neurons.




tFus: Transcranial Ultrasound Stimulation
rTMS : Repetitive Transcranial magnetic stimulation
tDCS : Transcranial Direct Current Stimulation (tDCS)
tACS : Transcranial Current Stimulation
tRNS : Transcranial Random Noise Stimulation

Image : <https://pubmed.ncbi.nlm.nih.gov/29311747/>

Repetitive Transcranial magnetic stimulation (rTMS)

- Produce changes in neuronal activity in regions of the brain.
- Help to detect critical region.
- Commonly used method.
- High temporal resolution.
- Help to detect particular behaviour by identifying brain region.
- Upper-limb motor-control recovery.
- Mood control and depression.
- Decreased activity in depression.



Approaches for EEG/MEG data analysis

- Time Domain Analysis : Event Related Potential
- Spectral analyses: amplitude, phase and coordination changes
- Event-related phenomena in comparison: phase-locked vs time-locked




Figure 9.1

Panel A shows single-trial EEG traces from 12 randomly selected trials (number inside plot indicates trial number). Data are from electrode FCz. Panel B shows 99 single trials in gray and their average—the ERP—in black. Panel C shows the same ERP with focused y-axis scaling.





Image :

<https://direct.mit.edu/books/book/4013/Analyzing-Neural-Time-Series-DataTheory-and, Ch 9>

Image : <https://www.nature.com/articles/nrn2774.pdf>

Transcranial Direct Current Stimulation (tDCS)

- Treat a range of neuropsychiatric disorders.
- Nature of stimulation either improve cognition or make it worse.
- Constant current not time varying
- Improved the clinical motor outcome of Parkinson's patients, as measured by the fine motor-control performance of their distal upper- limb.
- Used to treat patients with higher learning deficits.
- Enhanced language acquisition and fluency in healthy individuals.



Sensory Stimulation




Image: https://link.springer.com/chapter/10.1007/978-3-319-72038-8_16

- Multiple cortical areas specialised for sensory recognition or motor execution.
- Some functional modules are participating in both sensory as well as motor processing to analyze them sensory stimulation.
- Synchronized with EEG.

Transcranial Current Stimulation (tACS)

- Alternative Current
- Directly interfering with cortical rhythms (generating rhythms).
- Theta (improve cognition).
- Alpha (improve motor performance).
- Beta (deteriorate motor performance).
- Gamma (possibly interfere with attention).
 - Ability to detect motor behaviour, attention or visual stimuli
- Helps to diagnose Parkinson's disease or schizophrenia by attenuating or resetting anomalous oscillations.




Image: <https://www.mdpi.com/2076-3425/10/12/932>

tACS Continued..

- tACS stimulation improves visual performance when applied at theta-band over the occipital cortex, but not when stimulating over the medial prefrontal cortex or at alpha-band.
- Improved memory performance.
- Improved behavioural performance.
- Improved visual memory-matching reaction times.
- Modulated speech comprehension.
- Electrode positioning and stimulatory intensity are two important stimulatory parameters that impact the efficacy.

Adverse effects

- No severe negative effects have been reported in tDCS and tACS.
- TMS (seizures and auditory damage, ~5% incidence of mild adverse effects, including headache, nausea or local pain over the stimulated area).
- tDCS (Skin problems).
- tACS(none).

Conclusion

In sum, NIBS can modulate brain activity by enhancing the oscillatory coherence and synaptic plasticity, and, hence, has the potential to produce long-term therapeutic effects for different motor, language and memory disabilities.

Thank You!
Any Questions ??